


WinLIMS Data Download

You can access results from the same web-based application ([WinLIMS](#)) used to submit bottle orders or pre-log samples. Preliminary and approved results for all tests can be viewed using the 'Reports List' on the home page.



Clicking on the magnifying glass next to the [NELAC Level I.rpt](#) report will create a pop-up that will prompt you for your order number. After entering your order number (see above) you should see a PDF report similar to what the lab has generated in the past.

If you're looking for all **approved** data for your project you can click the magnifying glass next to [ApprovedExport – Project.rpt](#). The same pop-up window will prompt you for your project number and a date range. The generated report can be exported to a number of different formats including .csv and .xls by clicking the  button, selecting the desired file format and clicking [Export](#). A file named "CryViewer" will be generated and can usually be found in your downloads folder (or wherever your default download folder is).

To view **approved** data for all projects (ie: everything) use the [ApprovedExport – Program.rpt](#). The pop-up window will prompt you for your program number and a date range. Exporting the data is handled in the exact same manner as the Project report.