Common Laboratory Glassware

VT WSMD Wastewater Program Lab Manual Section #3

BURETS. Burets are used to deliver accurate volumes. They are especially useful in a procedure called titration.

FLASKS. Flasks are used for containing and mixing chemicals. There are







Support, Burette and Burette Clamp

Erlenmeyer Flask

Flat Bottom Flask







Burette

Boiling Round Bottom Flask

Filtering Flask









Automatic Burette

Volumetric Flask

Kjeldahl Flask

Distilling Flask

Common Laboratory Glassware

Bottles are used to store chemicals, to collect samples for testing purposes, and to dispense liquids.

Separatory funnels are used to separate one chemical mixture from another. The separated chemical usually is dissolved in one of two layers of liquid.







Reagent **Bottle**

BOD **Bottles**

Separatory **Funnel**

FUNNELS. A funnel is used for pouring solutions or transferring solid chemicals. This funnel also can be used with filter paper to remove solids from a solution.

TUBES. Test tubes are used for mixing small quantities of chemicals. They are also used as containers for bacterial testing (culture tubes).



Funnels

Test Tube



Culture Tube

Without Lip

A Buchner funnel is used to separate solids from a mixture. It is used with a filter flask and a vacuum.



Buchner Funnel



Wash Bottles

Common Laboratory Glassware

IMHOFF CONE. The Imhoff cone is used for the settleable solids testing of wastewater.




BEAKERS. Beakers are common pieces of laboratory equipment. They come in sizes from 1 ml to 4000 ml. They are used mainly for mixing chemicals and to measure approximate volumes.



GRADUATED CYLINDERS. Graduated cylinders

also are basic to any laboratory and come in sizes from 5 ml to 4000 ml. They are used to measure volumes more accurately than beakers.



Common Laboratory Glassware

OTHER LABWARE AND EQUIPMENT





Desiccator



Condenser







NIST Traceable Lab Thermometer







Volumetric Pipettes



Adjustable Volumetric Pipette

OTHER LABWARE AND EQUIPMENT (continued....)



Volumetric Flask Serological Pipettes